


The University of Iowa

At the Forefront of  
Free Radical Research

**Garry R. Buettner, Ph.D.**  
ESR Facility  
The University of Iowa  
Iowa City, IA 52242-1101  
[garry-buettner@uiowa.edu](mailto:garry-buettner@uiowa.edu)

**DETAPAC (DTPA), aqueous stock solution**

(1988, updated 2007-04-19; 2008-01-02; 2008.07.01)




**1. Introduction/Background**

Diethylenetriaminepentaacetic acid (DTPA or DETAPAC) is a metal chelating agent that was introduced to the free radical community in 1978 as a tool to modulate the catalytic activity redox active metals in superoxide-generating systems [1,2,]. It was later demonstrated that it is also effective in blunting the catalysis of ascorbate oxidation by both iron and copper [3,4,5,]. DETAPAC is a cousin of EDTA (See structures below.), but the Fe(III)DETAPAC complex is much more slowly reduced by superoxide (and ascorbate) than Fe(III)EDTA [6], **Table 1**.

**Table 1. Rate Constants for Reaction 1 [7,8,9]**

Ligand	$k_1/\text{M}^{-1} \text{ s}^{-1}$ (pH 7.4)
EDTA	$2 \times 10^6$
DETAPAC	$< 10^4$
HEDTA	$7.6 \times 10^5$
Desferal	$10^2\text{-}10^3$
Transferrin	$\leq 10^5$

DETAPAC (DTPA) is shown in the fully deprotonated form. However, it will not be fully deprotonated at near neutral pH values. The usual suppliers sell DTPA as the penta-acid. As such it is very difficult to dissolve and produces a very acidic solution. For most biochemical and biological experiments we add some NaOH to the stock solution to partially neutralize the DTPA, producing a near neutral solution that is much more easily dissolved and will not challenge the pH of the solution to which it is being added.



## 2. Materials

Source: Example, Aldrich 09,390-2      Sigma D-1133 (1 kg )

FW = 393.35 Da for the pentaacid,  $C_{14}H_{23}N_3O_{10}$  #CAS 67-43-6

The acid dissociation constants for DTPA are [10, 11]:

H <sub>5</sub> A	pK <sub>1</sub> = 1.79
	pK <sub>2</sub> = 2.56
	pK <sub>3</sub> = 4.42
	pK <sub>4</sub> = 8.76
	pK <sub>5</sub> = 10.42

## 3. Procedure:

To make 50 mL of 10 mM stock solution of DETAPAC:

$$\frac{393.35 \text{ g}}{\text{mole}} \times \frac{0.01 \text{ mole}}{1 \text{ L}} \times 0.05 \text{ L} \times \frac{1000 \text{ mg}}{1 \text{ g}} = 196.67 \text{ mg}$$

How much NaOH should we add to dissolve the DTPA easily and get a near-neutral solution?

Looking at the pK<sub>a</sub>s it would appear we should add about 3.4 equivalents of base per mole of DETAPAC. In this stock solution (0.01 M, 50 mL) we have a total of 0.5 mmole of DETPAC. Thus, we should add:

0.5 mmol × 3.4 = 17 mmol of NaOH; this much NaOH would be contained in 1.7 mL of 1 M NaOH, or 17 mL of 0.1 M NaOH.

Thus, to make 50 mL of a 0.01 M solution of DETAPAC

1. To a 50 mL volumetric flask add 1.7 mL of 1.0 M NaOH solution (17 mL of 0.1 M);
2. add 197 mg of DETAPAC, pentaacid form;
3. add clean water (DI, Nanopure, Milli-Q, whatever) for a total volume of approximately 40-45 mL;
4. dissolve by stirring or sonication, we prefer sonication;
5. adjust pH to approximately 7, or 7.4 or whatever is best for most proposed uses;
6. dilute to mark.

Scale up or down, of course, for other volumes or concentrations for stock solutions. DETAPAC is typically used at 0.1 mM or less in experiments, so a 10 mM stock is convenient.

**4. Hazards:** The bottle of solid DTPA and MSDS for DTPA have many hazardous warnings. However, EDTA comes with these same warnings. We eat EDTA and DTPA is given *i.v.* to remove heavy metals, e.g. plutonium. Thus, prudent handling is in order, but perspective should be maintained.

## 5. References:

---

- 1 Buettner GR, Oberley LW, Leuthauser SWHC. (1978) The effect of iron on the distribution of superoxide and hydroxyl radicals as seen by spin trapping and on the superoxide dismutase assay." *Photochem Photobiol* **28**: 693-695.  
<http://www.healthcare.uiowa.edu/CoreFacilities/esr/publications/buettnerpubs/pdf/P&P-1978-28-693-GRB-Fe-SOD.pdf>
- 2 Buettner GR, Oberley LW. (1978) Considerations in the spin trapping of superoxide and hydroxyl radicals in aqueous systems using 5,5-dimethyl-1-pyrroline-1-oxide. *Biochem Biophys Res Commun* **83**: 69-74. [<http://www.healthcare.uiowa.edu/CoreFacilities/esr/contact/buettner/pdf/BBRC-1978-83-69-DMPO.pdf>] or [[http://dx.doi.org/10.1016/0006-291X\(78\)90398-4](http://dx.doi.org/10.1016/0006-291X(78)90398-4)]
- 3 Buettner GR. (1986) Ascorbate autoxidation in the presence of iron and copper chelates. *Free Rad Res Commun* **1**: 349-353.  
<http://www.healthcare.uiowa.edu/CoreFacilities/esr/publications/buettnerpubs/pdf/FRR-1986-1-349-grb-Asch-metals.pdf>
- 4 Buettner GR. (1988) In the absence of catalytic metals, ascorbate does not autoxidize at pH 7: Ascorbate as a test for catalytic metals. *J Biochem Biophys Meth* **16**: 20-40.  
[http://www.healthcare.uiowa.edu/corefacilities/esr/publications/buettnerpubs/pdf/JBBM-1988-16-27-GRB-Ascorbate\\_test.pdf](http://www.healthcare.uiowa.edu/corefacilities/esr/publications/buettnerpubs/pdf/JBBM-1988-16-27-GRB-Ascorbate_test.pdf) or [[http://dx.doi.org/10.1016/0165-022X\(88\)90100-5](http://dx.doi.org/10.1016/0165-022X(88)90100-5)]
- 5 Buettner GR, Jurkiewicz BA. (1996) Catalytic metals, ascorbate, and free radicals: combinations to avoid. *Radiation Research* **145**:532-541.  
<http://www.healthcare.uiowa.edu/corefacilities/esr/publications/buettnerpubs/pdf/RadRes-1996-145-532-Buettner-Acs.pdf> or <http://www.jstor.org/stable/3579271>
- 6 Buettner GR, Doherty TP, Patterson LK. (1983) The kinetics of the reaction of superoxide with Fe(III)EDTA, Fe(III)DETAPAC and Fe(III)HEDTA. *FEBS Lett* **183**: 143-146.  
[<http://www.healthcare.uiowa.edu/corefacilities/esr/publications/buettnerpubs/pdf/13FEBS-1983-158-143-FeDETAPAC.pdf>] or [[http://dx.doi.org/10.1016/0014-5793\(83\)80695-4](http://dx.doi.org/10.1016/0014-5793(83)80695-4)]
- 7 Ross A.B., Mallard W.G., Helman W.P., Buxton G.V., Huie R.E., Neta P. (1994) NDRL-NIST Solution Kinetics Database: Ver. 2.0. National Institute of Standards and Technology, Gaithersburg. The modern version can be accessed via the web at;  
<http://www.rad.nd.edu/rcdc/index.html>
- 8 Buettner G.R., Doherty T.P., Patterson L.K. (1983) The kinetics of the reaction of superoxide with Fe(III)DETAPAC and Fe(III)HEDTA. *FEBS Lett.* **183**:143-146.  
[<http://www.healthcare.uiowa.edu/corefacilities/esr/publications/buettnerpubs/pdf/13FEBS-1983-158-143-FeDETAPAC.pdf>] or [[http://dx.doi.org/10.1016/0014-5793\(83\)80695-4](http://dx.doi.org/10.1016/0014-5793(83)80695-4)]
- 9 Buettner, G.R. (1987) Activation of oxygen by metal complexes and its relevance to autoxidative processes in living systems. *Bioelectrochem. Bioenergetics*, **18**: 29-36.  
[<http://www.healthcare.uiowa.edu/corefacilities/esr/publications/buettnerpubs/pdf/BioelectrochemBioen-1987-18-29-grb-Metals.pdf>] or [[http://dx.doi.org/doi:10.1016/0302-4598\(87\)85005-5](http://dx.doi.org/doi:10.1016/0302-4598(87)85005-5)]
- 10 Sille'n LG, Martell AE. (1964) *Stability Constants of Metal-Ion Complexes*, Special Publication No. 17; The Chemical Society: London.
- 11 Sille'n LG, Martell AE. (1971) *Stability Constants of Metal-Ion Complexes Supplement No. 1*, Special Publicatcon No. 25; The Chemical Society: London.