



Product application guide

Prepare the area: Cleanse intact skin. Dry the surrounding skin thoroughly. Ensure that skin is free of dimethicone, skin sealants and emollients. Use of skin barrier under dressing is not necessary.




1. Area to protect. Assess the patient's anatomy and determine appropriate dressing positioning.



2. After skin is prepared, remove the center release film by gently pulling on pink-lined edge.



3. Hold buttocks apart. Apply dressing to sacral area and into upper aspect of gluteal cleft, with dressing "base" positioned to cover coccyx area.



4. Remove side release films and gently smooth each side into place.



5. Product placement.



6. Press and smooth the dressing to ensure the entire dressing is in contact with the skin.



Proper fit:

- Dressing should be positioned to fit patient anatomy.
- Ensure dressing conforms to the skin and avoid gaps or air pockets.

Pressure ulcer prevention re-application guide



1. Assess to confirm dressing is intact and applied correctly.



2. Gently pull handling tabs to begin to release dressing from skin.



3. Continue to release dressing from skin using handling tabs until skin exposed for skin check.



4. While maintaining dressing position at gluteal cleft, perform assessment of skin.



5. Reapply the foam and borders of the dressing.



6. Confirm dressing is replaced to its original position, making sure the border is intact and flat.



7. Press and smooth the dressing to ensure the entire dressing is in contact with the skin.